AYUMiEYE

medical

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Preprints

HOW IT WORKS
INSTRUCTIONS FOR AUTHORS

SUBJECT AREAS
ADVISORY BOARD

SCREENING PREPRINTS
ABOUT

preprints.org > doi: 10.20944/preprints202012.0336.v1

Preprint Article Version 1 Preserved in Portico This version is not peer-reviewed

Comparison of Gait Analysis Between a Triaxial Accelerometer-Based Device and an Optical Motion Capture System

Yutaka Ito * and Yuichi Ota

Version 1 : Received: 11 December 2020 / Approved: 14 December 2020 / Online: 14 December 2020 (13:22:52 CFT)

How to cite: Ito, T.; Ota, Y. Comparison of Gait Analysis Between a Triaxial Accelerometer-Based Device and an Optical Motion Capture System. *Preprints* 2020, 2020120336 (doi: 10.20944/preprints202012.0336.v1).

Abstract

AYUMI EYE is an accelerometer-based gait analysis device that measures the 3D accelerations of the human trunk. This study investigated the measurement accuracy of the AYUMI EYE as hardware as well as the accuracy of the gait cycle extraction program via simultaneous measurements using AYUMI EYE, a ground reaction force (GRF), and an optical motion capture system called VICON. The study was conducted with four healthy individuals as participants. The gait data were obtained by simulating four different patterns for three trials each: normal walking, anterior-tilt walking, hemiplegic walking, and shuffling walking. The AYUMI EYE and VICON showed good agreement for both the acceleration and displacement data. The durations of subsequent stride cycles calculated using the AYUMI EYE and GRF were in good agreement based on the calculated cross-correlation coefficients (CCs) with an r value of 0.896 and p-value less than 0.05, and their accuracies for these results were sufficient.

Subject Areas

AYUMI EYE; acceleration sensor; gait analysis; optical sensor

Copyright: This is an open access article distributed under the [Creative Commons Attribution License](#) which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.

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